

EXPANDERTM

ESOPHAGEAL DILATOR SET



EXPANDERTM

ESOPHAGEAL DILATOR SET

Product Features : -

- Set of Six & Eight Dilators
- Total Length 85 cm
- Marking Every 1 cm
- Built - in Radiopaque Metal Ring
- Highly Smooth and Flexible

Catalog No.	Set Includes	Dilator Sizes
EDS - 0106085	Six Dilators + One Spring Tip Wire	5,7,9,11,13,15,mm x 85 cm
EDS - 0108065	Eight Dilators + One Spring Tip Wire	5,7,9,11,13,15,17,19 mm x 85 cm

Other size and configurations are available on request.

Manufactured by:

 **Endo-Med[®]**
Technologies Pvt.Ltd.

New Delhi - INDIA
Ph. : +91 11 22048222
Fax : +91 11 22447748
URL : www.endomed.in
E-mail : info@endomed.in

Distributed by: