

PALLIATORTM

ESOPHAGEAL STENT


PALLIATOR



PALLIATOR
(Anti-Reflux)



Scaled making on delivery System
helps deployment without floro.




PALLIATORTM

ESOPHAGEAL STENT

Product Features : -

- Nitinol Stent Highly Radiopaque.
- Covered From Inside.
- Easy Retrieval using Drawstring.
- Dynamic Radial Force for Maximum Lumen Size
- Dumble Type for Anti-Migration
- Flexible Construction to Conform to Anatomy
- X-ray Visible Marks
- Scaled Marking on Delivery System



Catalog No.	Stent OD	Stent Length	Catheter Size
EMS - 0118008	18 mm	8 cm	7 mm x 80 cm
EMS - 0118010	18 mm	10 cm	7 mm x 80 cm
EMS - 0118012	18 mm	12 cm	7 mm x 80 cm
EMS - 0118015	18 mm	15 cm	7 mm x 80 cm
EMS - 0118010 (AR)	18 mm	10 cm	7 mm x 80 cm
EMS - 0118012 (AR)	18 mm	12 cm	7 mm x 80 cm
EMS - 0118015 (AR)	18 mm	15 cm	7 mm x 80 cm

Other size and configurations are available on request.

Manufactured by:

 **Endo-Med[®]**
Technologies Pvt.Ltd.

New Delhi - INDIA
Ph. : +91 11 22048222
Fax : +91 11 22447748
URL : www.endomed.in
E-mail : info@endomed.in

Distributed by: